



The Difference

- Certified to TOC Levels < 30ppb
- Pipettor like disconnect to eliminate cutting
- Excellent solvent compatibility for HPLC & IMS testing
- Designed to minimize risks of cross contamination

Physical Properties

Part	Material
Handle	Polypropylene
Tip	Knitted Polyester
Tip Shaft	Polypropylene
Tip Tray	PET

The SurfCHECK® Cleaning Validation Swab is a unique, patent-pending tool specifically engineered to support pharmaceutical cleaning validation work using various analytical methods, including TOC, HPLC, IMS & UV-Vis. The SurfCHECK® will eliminate background interference from improper technician handling techniques by a design that allows technicians to use the tools the same way every time without ever touching the swab tip.

Each tip is firmly held in place in an easy-to-use tray, connecting to the handle at point-of-use. Tips are oriented for ease of use, ensuring that technicians will insert them properly and securely. After swabbing, the tips can be placed into the diluent by depressing a protected button at the top of the pen. There is no need for cutting implements or snapping mechanisms which can contaminate the sample. Each tip and SurfCHECK® are certified at less than 30ppb TOC for that type of analytical work.

Scan Here



The Next Generation Cleaning Validation Swab

- 1** Hands Free Attachable Disposable Heads
- 2** Extremely Flexible Tip
- 3** Detaches with a Click of a Button
- 4** Ready for validation, Prevents Cross Contamination

PART NO.	DESCRIPTION	UNITS/CASE
HT200A	Flexible Medium Cylinder Polyester Tip Swabs, 2 FTSWAB-PENS, 6 inch Handle - 12 Swabs/Tray, 12 Trays/Kit	12 Kits



**Foamtec
International**

Wilshire Contamination Control Division

1.866.922.7160 (USA Toll-Free)
 1.760.599.7160 (USA & Europe)
 66.35740717 (Thailand)
www.foamtecintlwcc.com